

Product datasheet for RC212216

PRAS40 (AKT1S1) (NM_001098633) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PRAS40 (AKT1S1) (NM_001098633) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PRAS40
Synonyms:	Lobe; PRAS40
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC212216 representing NM_001098633 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGTCGGGGCGCCCCGAGGAGCTGTGGGAGGCCGTGGTGGGGCCGCTGAGCGCTTCCGGGCCGGA
CTGGCACGGAGCTGGTGTGCTGACCGCGGCCCGCCGCCACCACCCGCCCCGGGCCCTGTGCCTATGC
TGCCCATGGTCGAGGAGCCCTGGCGGAGGCAGCGCCGTTGCCTCCACGACATCGCACTGGCCACAGG
GCTGCCACTGCTGCTCGGCTCCTGCGCCCCACCAGCACACAGCCACCCAGTCCACACCCAGCCAC
CCCGGCCTACCTGGCCAGAGAGACAACGAGGAGGACGAGGATGAGCCACAGACAGAGACCTCCGG
GGAGCAGCTGGGCATTAGTGATAATGGAGGGCTCTTTGTGATGGATGAGGACGCCACCTCCAGGACCTT
CCCCCTTCTGTGAGTCAGACCCGAGAGTACAGATGATGGCAGCCTGAGCGAGGAGACCCCGCGGCC
CCCCACCTGCTCAGTGCCCCAGCCTCAGCCCTACCCACACAGCAGTACGCCAAGTCCCTGCCTGTGTC
TGTGCCCGTCTGGGGCTTCAAGGAGAAGAGGACAGAGGCGCGGTTCATCAGATGAGGAGAAATGGGCCGCC
TCTTCGCCGACCTGGACCGCATCGCGGCAGCATGCGCGCGTGGTGTGCGAGAGGCCGAGGACACCC
AGGTCTTCGGGGACCTGCCACGGCCGCGGCTTAACACCAGCGACTTCCAGAAGCTGAAGCGGAAATAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC212216 representing NM_001098633
Red=Cloning site Green=Tags(s)

MASGRPEELWEAVVGAERFRARTGTELVLLTAAPPPPPRPGPCAYAAHGRGALAEAARRCLHDIALAHR
 AATAARPPAPPAPQPPSPTSPPPRPTLAREDNEEDEDEPTETETETSGEQLGSDNGLFVMDATLQDL
 PPFCESDPESTDDGSLSEETPAGPPTCSVPPASALPTQQYAKSLPVSVPVWGFKEKRTEARSSDEENGPP
 SSPDLDRIAASMRALVLEAEDTQVFGDLPRPRLNTSDFQKLKRKY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8025_h05.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001098633

ORF Size: 768 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001098633.4](#)

RefSeq Size: 1789 bp

RefSeq ORF: 771 bp

Locus ID: 84335


UniProt ID: [Q96B36](#)

Cytogenetics: 19q13.33


MW: 27.2 kDa

Gene Summary: AKT1S1 is a proline-rich substrate of AKT (MIM 164730) that binds 14-3-3 protein (see YWHAH, MIM 113508) when phosphorylated (Kovacina et al., 2003 [PubMed 12524439]).[supplied by OMIM, Mar 2008]

Product images:



Circular map for RC212216



Coomassie blue staining of purified AKT1S1 protein (Cat# [TP312216]). The protein was produced from HEK293T cells transfected with AKT1S1 cDNA clone (Cat# RC212216) using MegaTran 2.0 (Cat# [TT210002]).