

## Product datasheet for **RC212222L4V**

### HEPACAM2 (NM\_198151) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	HEPACAM2 (NM_198151) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HEPACAM2
Synonyms:	MIKI
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_198151
ORF Size:	1350 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212222).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_198151.1</a>
RefSeq Size:	2110 bp
RefSeq ORF:	1353 bp
Locus ID:	253012
UniProt ID:	<a href="#">A8MWW5</a>
Cytogenetics:	7q21.2
Protein Families:	Transmembrane
MW:	49.9 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein related to the immunoglobulin superfamily that plays a role in mitosis. Knockdown of this gene results in prometaphase arrest, abnormal nuclear morphology and apoptosis. Poly(ADP-ribosylation) of the encoded protein promotes its translocation to centrosomes, which may stimulate centrosome maturation. A chromosomal deletion including this gene may be associated with myeloid leukemia and myelodysplastic syndrome in human patients. [provided by RefSeq, Oct 2016]