

Product datasheet for RC212227

OR4C15 (NM_001001920) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OR4C15 (NM_001001920) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	OR4C15
Synonyms:	OR11-127; OR11-134
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC212227 representing NM_001001920 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTTCTCAATGACAACAGAAGCACTCAATAATTTGCACTTGGATGTACCAACTGTTAATGACTATGA
TACCACAAATTGATCTGAAGCAAATTTCTTTGCTAATTGCAGACTATACATGATCCCTGTTGGAGC
TTTCATCTTTTCCTGGGAAACATGCAAACCAAAGCTTGTAACTGAGTTTGCCTCCTGGGACTTCA
CAGAAATCAAATGTTTCAGGAAATAGTATTTGTTGTATTTTGTGTTCTACATTGCAACTGTTGGGGCA
ACATGCTAATTGTAGTAACCATTCTCAGCAGCCCTGCTCTCTGGTGTCTCCTATGACTTCTTCTGGG
CTTCTGTCTTCTGGATGCGTGCTTCTCATCTGTCATCACCCAAAGATGATTGTAGACTCCCTCTAT
GTGACAAAAACCATCTCTTTGAAGGCTGCATGATGCAGCTCTTTGCTGAACACTTCTTGTGGGGTGG
AGGTGATTGTCCTCACAGCCATGGCCTATGATCGTTATGTGGCCATTTGCAAGCCCTTGCACTTCTTC
TATCATGAACAGGAGGCTCTGTGGCATTCTGATGGGGTAGCCTGGACAGGGGGCCTCTTGCACTCCATG
ATACAAATCTTTTTACTTTCCAGCTCCCTTTGTGGCCCAATGTCATCAATCACTTTATGTGTGACT
TGTACCGTTACTGGAGCTTGCCTGCAGTACTCACATCTTTGGCCTCATGGTGGTCAACAGTGG
GTTTATCTGCATCATAAATTCTCCTTGTGCTGTCTCCTATGCTGTCATCTTGTCTCTGAGAAC
CACAGTCTGAAGGGCGCTGAAAGCTCTCCACCTGTGGATCTCACATTGCTGTTGATTTTGTCT
TTGTCCCATGCATATTTGTATATACAGCCTCCATCTGCTTTTCCCTTGACAAAATGGCGGCAATATT
TTATATCATCTTAAATCCCTTGCTCAATCCTTTGATTTACACTTTCAGGAATAAGGAAGTAAACAGGCC
ATGAGGAGAATATGGAACAGACTGATGGTGGTTTCTGATGAGAAAAGAAAATATAAACTT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTAA



[View online »](#)

Protein Sequence: >RC212227 representing NM_001001920
 Red=Cloning site Green=Tags(s)

MFSMTTEALNNFALGCTNLLMTMIPQIDLKQIFLCPNCRLYMIPVGAFIFSLGNMQNSFVTEFVLLGLS
 QNPNVQEI V F V V F L F V Y I A T V G G N M L I V V T I L S S P A L L V S P M Y F L G F L S F L D A C F S S V I T P K M I V D S L Y
 V T K T I S F E G C M M Q L F A E H F F A G V E V I V L T A M A Y D R Y V A I C K P L H Y S S I M N R R L C G I L M G V A W T G L L H S M
 I Q I L F T F Q L P F C G P N V I N H F M C D L Y P L L E L A C T D T H I F G L M V V I N S G F I C I I N F S L L L V S Y A V I L L S L R T
 H S S E G R W K A L S T C G S H I A V V I L F F V P C I F V Y T R P P S A F S L D K M A A I F Y I I L N P L L N P L I Y T F R N K E Y K Q A
 M R R I W N R L M V V S D E K E N I K L

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms: https://cdn.origene.com/chromatograms/mk8001_d01.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001001920
ORF Size:	1110 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001001920.1 , NP_001001920.1
RefSeq Size:	1113 bp
RefSeq ORF:	951 bp
Locus ID:	81309
UniProt ID:	Q8NGM1
Cytogenetics:	11q11
Protein Families:	Transmembrane
Protein Pathways:	Olfactory transduction
MW:	41.5 kDa
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]