

Product datasheet for **RC212231L4V**

SMARCA6 (HELLS) (NM_018063) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SMARCA6 (HELLS) (NM_018063) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SMARCA6
Synonyms:	ICF4; LSH; Nbla10143; PASG; SMARCA6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018063
ORF Size:	2514 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212231).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018063.3
RefSeq Size:	3165 bp
RefSeq ORF:	2517 bp
Locus ID:	3070
UniProt ID:	Q9NRZ9
Cytogenetics:	10q23.33
Domains:	SNF2_N, DEAD, helicase_C
MW:	96.9 kDa



[View online »](#)

Gene Summary:

This gene encodes a lymphoid-specific helicase. Other helicases function in processes involving DNA strand separation, including replication, repair, recombination, and transcription. This protein is thought to be involved with cellular proliferation and may play a role in leukemogenesis. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jan 2014]