

## Product datasheet for RC212237L2V

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# NDRG2 (NM 201537) Human Tagged ORF Clone Lentiviral Particle

#### **Product data:**

**Product Type: Lentiviral Particles** 

**Product Name:** NDRG2 (NM 201537) Human Tagged ORF Clone Lentiviral Particle

Symbol: SYLD Synonyms:

**Mammalian Cell** None

Selection:

Vector: pLenti-C-mGFP (PS100071)

mGFP Tag:

ACCN: NM 201537 **ORF Size:** 1113 bp

**ORF Nucleotide** 

The ORF insert of this clone is exactly the same as(RC212237).

Sequence: OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeq: NM 201537.1, NP 963831.1

RefSeq Size: 2075 bp RefSeq ORF: 1116 bp Locus ID: 57447 **UniProt ID: Q9UN36** Cytogenetics: 14q11.2

MW: 40.6 kDa







### **Gene Summary:**

This gene is a member of the N-myc downregulated gene family which belongs to the alpha/beta hydrolase superfamily. The protein encoded by this gene is a cytoplasmic protein that may play a role in neurite outgrowth. This gene may be involved in glioblastoma carcinogenesis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2017]