

Product datasheet for **RC212255L3V**

NADE (BEX3) (NM_014380) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NADE (BEX3) (NM_014380) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NADE
Synonyms:	Bex; DXS6984E; HGR74; NADE; NGFRAP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_014380
ORF Size:	333 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212255).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014380.1
RefSeq Size:	984 bp
RefSeq ORF:	336 bp
Locus ID:	27018
UniProt ID:	Q00994
Cytogenetics:	Xq22.2
Domains:	BEX
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Neurotrophin signaling pathway

MW: 13 kDa

Gene Summary: May be a signaling adapter molecule involved in p75NTR-mediated apoptosis induced by NGF. Plays a role in zinc-triggered neuronal death (By similarity). May play an important role in the pathogenesis of neurogenetic diseases.[UniProtKB/Swiss-Prot Function]