

Product datasheet for **RC212263L3V**

DPS1 (PDSS1) (NM_014317) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DPS1 (PDSS1) (NM_014317) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DPS1
Synonyms:	COQ1; COQ1A; COQ10D2; DPS; hDPS1; SPS; TPRT; TPT; TPT 1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_014317
ORF Size:	1245 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212263).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014317.3
RefSeq Size:	1679 bp
RefSeq ORF:	1248 bp
Locus ID:	23590
UniProt ID:	Q5T2R2
Cytogenetics:	10p12.1
Protein Pathways:	Terpenoid backbone biosynthesis
MW:	46.1 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene is an enzyme that elongates the prenyl side-chain of coenzyme Q, or ubiquinone, one of the key elements in the respiratory chain. The gene product catalyzes the formation of all trans-polyprenyl pyrophosphates from isopentyl diphosphate in the assembly of polyisoprenoid side chains, the first step in coenzyme Q biosynthesis. The protein may be peripherally associated with the inner mitochondrial membrane, though no transit peptide has been definitively identified to date. Defects in this gene are a cause of coenzyme Q10 deficiency. [provided by RefSeq, Jul 2008]