

Product datasheet for **RC212310A1V**

Human LILRA5 (NM_021250) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human LILRA5 (NM_021250) AAV Particle
Tag:	Myc-DDK
Symbol:	LILRA5
Synonyms:	CD85; CD85F; ILT-11; ILT11; LILRB7; LIR-9; LIR9
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC212310 representing NM_021250 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCACCATGGTCTCATCCATCTGCACAGCTGCAGCCAGTGGGAGGAGACGCCGTGAGCCCTGCCCTCA
TGGTCTGCTCTGCCTCGGGCTGAGTCTGGGCCCCAGGACCCACGTGCAGGCAGGGAACCTCTCCAAAGC
CACCTCTGGGCTGAGCCAGGCTCTGTGATCAGCCGGGGAACTCTGTGACCATCCGGTGTGAGGGGACC
CTGGAGGCCAGGAATACCGTCTGGTTAAAGAGGGAAGCCAGAACCTGGGACACAGAACCCACTGG
AGCCCAAGAACAAGGCCAGATTCTCCATCCATCCATGACAGAGCACCATGCAGGGAGATACCGCTGTTA
CTACTACAGCCCTGCAGGCTGGTCAGAGCCAGCGACCCCTGGAGCTGGTGGTGACAGGATTCTACAAC
AAACCCACCCTCTCAGCCCTGCCAGTCTGTGGTGACCTCAGGAGAGAACGTGACCCCTCAGTGTGGCT
CACGGCTGAGATTCGACAGGTTCACTCTGACTGAGGAAGGAGACCACAAGCTCTCCTGGACCTTGGACTC
ACAGCTGACCCCAAGTGGGCACTTCCAGGCCCTGTTCCCTGTGGGCCCTGTGACCCCAAGCCACAGGTGG
ATGCTCAGATGCTATGGCTCTCGCAGGCATATCCTGCAGGTATGGTCAGAACCCAGTACCTCTGGAGA
TTCCGGTCTCAGGAGCAGCTGATAACCTCAGTCCGTACAAAACAAGTCTGACTCTGGGACTGCCTCACA
CCTTCAGGATTACGCAGTAGAGAATCTCATCCGATGGGCATGGCCGGCTTGATCTGGTGGTCTTGGG
ATTCTGATATTTAGGATTGGCACAGCCAGAGAAGCCCCAAGCTGCAGCTGGAAGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC212310 representing NM_021250 Red=Cloning site Green=Tags(s) MAPWSPSAQLQPVGGDAVSPALMVLLCLGLSLGPRTHVQAGNLSKATLWAEPGSVISRGNSVTIRCQGT LEAQEYRLVKEGSPEPWDTONPLEPKNKARFSIPSMTEHHAGRYRCYYSPAGWSEPSDPLELVVTGFYN KPTLSALPSPVVTSGENVTLQCGSRLRFDRFILTEEGDHKLSWTLDSQLTPSGQFQALFPVGPVTPSHRW MLRCYGSRRHILQVWSEPSDLLEIPVSGAADNLSPSQNKSDSGTASHLQDYAVENLIRMGAGLILVVLG ILIFQDWHSQRSPQAAAGR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_021250
ORF Size:	897 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_021250.2
RefSeq Size:	1365 bp
RefSeq ORF:	900 bp
Locus ID:	353514
UniProt ID:	A6NI73
Cytogenetics:	19q13.42
MW:	32.6 kDa