

Product datasheet for RC212315L4

ENTPD7 (NM_020354) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ENTPD7 (NM_020354) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	ENTPD7
Synonyms:	LALP1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212315).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_020354
ORF Size:	1812 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_020354.2
RefSeq Size:	8564 bp
RefSeq ORF:	1815 bp
Locus ID:	57089
UniProt ID:	Q9NQZ7
Cytogenetics:	10q24.2
Domains:	GDA1_CD39
Protein Families:	Transmembrane
MW:	69 kDa
Gene Summary:	This gene encodes a purine-converting ectoenzyme which belongs to the ecto-nucleoside triphosphate diphosphohydrolase (E-NTPDase) family. The encoded protein hydrolyzes extracellular nucleoside triphosphates (UTP, GTP, and CTP) to nucleoside monophosphates as part of a purinergic signaling pathway. It contains two transmembrane domains at the N- and C-termini and a large, hydrophobic catalytic domain located in between. This gene affects oxidative stress as well as DNA damage and is a mediator of senescence. [provided by RefSeq, Mar 2017]

Product images:



Circular map for RC212315L4