

Product datasheet for **RC212317L3V**

PPCS (NM_024664) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PPCS (NM_024664) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PPCS
Synonyms:	CMD2C
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_024664
ORF Size:	933 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212317).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024664.2
RefSeq Size:	2312 bp
RefSeq ORF:	936 bp
Locus ID:	79717
UniProt ID:	Q9HAB8
Cytogenetics:	1p34.2
Domains:	dfp
Protein Pathways:	Metabolic pathways, Pantothenate and CoA biosynthesis



[View online »](#)

MW: 34 kDa

Gene Summary: Biosynthesis of coenzyme A (CoA) from pantothenic acid (vitamin B5) is an essential universal pathway in prokaryotes and eukaryotes. PPCS (EC 6.3.2.5), one of the last enzymes in this pathway, converts phosphopantothenate to phosphopantothenoylcysteine (Daugherty et al., 2002 [PubMed 11923312]).[supplied by OMIM, Mar 2008]