

## **Product datasheet for RC212340L1**

#### 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com

OriGene Technologies, Inc.

techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### ATG16L1 (NM\_017974) Human Tagged Lenti ORF Clone

**Product data:** 

**Product Type:** Expression Plasmids

Product Name: ATG16L1 (NM 017974) Human Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: ATG16L1

Synonyms: APG16L; ATG16A; ATG16L; IBD10; WDR30

Mammalian Cell None

Selection:

Vector:pLenti-C-Myc-DDK (PS100064)E. coli Selection:Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC212340).

Sequence:

**Restriction Sites:** Sgfl-Mlul

**Cloning Scheme:** 





<sup>\*</sup> The last codon before the Stop codon of the ORF

**ACCN:** NM\_017974

ORF Size: 1764 bp





#### ATG16L1 (NM\_017974) Human Tagged Lenti ORF Clone - RC212340L1

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

**RefSeq:** <u>NM 017974.3</u>

 RefSeq Size:
 3354 bp

 RefSeq ORF:
 1767 bp

 Locus ID:
 55054

 UniProt ID:
 Q676U5

 Cytogenetics:
 2q37.1

Domains: WD40

MW: 66.1 kDa


**Gene Summary:** The protein encoded by this gene is part of a large protein complex that is necessary for

autophagy, the major process by which intracellular components are targeted to lysosomes for degradation. Defects in this gene are a cause of susceptibility to inflammatory bowel disease type 10 (IBD10). Several transcript variants encoding different isoforms have been


found for this gene.[provided by RefSeq, Jun 2010]



# **Product images:**



Circular map for RC212340L1



Double digestion of RC212340L1 using Sgfl and Mlul  $\,$