

## Product datasheet for **RC212340L3V**

### **ATG16L1 (NM\_017974) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	ATG16L1 (NM_017974) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ATG16L1
Synonyms:	APG16L; ATG16A; ATG16L; IBD10; WDR30
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017974
ORF Size:	1764 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212340).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_017974.3</a>
RefSeq Size:	3354 bp
RefSeq ORF:	1767 bp
Locus ID:	55054
UniProt ID:	<a href="#">Q676U5</a>
Cytogenetics:	2q37.1
Domains:	WD40
MW:	66.1 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is part of a large protein complex that is necessary for autophagy, the major process by which intracellular components are targeted to lysosomes for degradation. Defects in this gene are a cause of susceptibility to inflammatory bowel disease type 10 (IBD10). Several transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jun 2010]