

Product datasheet for **RC212368L2V**

Menin (MEN1) (NM_130799) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Menin (MEN1) (NM_130799) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Menin
Synonyms:	MEAI; SCG2
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_130799
ORF Size:	1830 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212368).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_130799.1
RefSeq Size:	2772 bp
RefSeq ORF:	1833 bp
Locus ID:	4221
UniProt ID:	O00255
Cytogenetics:	11q13.1
Domains:	Menin
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)

MW: 67.3 kDa

Gene Summary: This gene encodes menin, a tumor suppressor associated with a syndrome known as multiple endocrine neoplasia type 1. Menin is a scaffold protein that functions in histone modification and epigenetic gene regulation. It is thought to regulate several pathways and processes by altering chromatin structure through the modification of histones. [provided by RefSeq, May 2019]