



<b>Serotype:</b>	AAV-2
<b>ACCN:</b>	NM_182503
<b>ORF Size:</b>	573 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_182503.2</a></u>
<b>RefSeq Size:</b>	6270 bp
<b>RefSeq ORF:</b>	576 bp
<b>Locus ID:</b>	134637
<b>UniProt ID:</b>	<u><a href="#">Q7Z6V5</a></u>
<b>Cytogenetics:</b>	6q24.2
<b>MW:</b>	20.9 kDa