

## Product datasheet for **RC212397L2V**

### IL20 (NM\_018724) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	IL20 (NM_018724) Human Tagged ORF Clone Lentiviral Particle
Symbol:	IL20
Synonyms:	IL-20; IL10D; ZCYTO10
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_018724
ORF Size:	528 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212397).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_018724.3</a> , <a href="#">NP_061194.2</a>
RefSeq Size:	1252 bp
RefSeq ORF:	531 bp
Locus ID:	50604
UniProt ID:	<a href="#">Q9NYY1</a>
Cytogenetics:	1q32.1
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway



[View online »](#)

**MW:** 20.1 kDa

**Gene Summary:** The protein encoded by this gene is a cytokine structurally related to interleukin 10 (IL10). This cytokine has been shown to transduce its signal through signal transducer and activator of transcription 3 (STAT3) in keratinocytes. A specific receptor for this cytokine is found to be expressed in skin and upregulated dramatically in psoriatic skin, suggesting a role for this protein in epidermal function and psoriasis. [provided by RefSeq, Jul 2008]