

Product datasheet for **RC212434L2V**

ASAH1 (NM_004315) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ASAH1 (NM_004315) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ASAH1
Synonyms:	AC; ACDase; ASAH; PHP; PHP32; SMAPME
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_004315
ORF Size:	1233 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212434).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004315.2
RefSeq Size:	2503 bp
RefSeq ORF:	1236 bp
Locus ID:	427
UniProt ID:	Q13510
Cytogenetics:	8p22
Domains:	CBAH
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Lysosome, Metabolic pathways, Sphingolipid metabolism

MW: 46.3 kDa

Gene Summary: This gene encodes a member of the acid ceramidase family of proteins. Alternative splicing results in multiple transcript variants, at least one of which encodes a preproprotein that is proteolytically processed. Processing of this preproprotein generates alpha and beta subunits that heterodimerize to form the mature lysosomal enzyme, which catalyzes the degradation of ceramide into sphingosine and free fatty acid. This enzyme is overexpressed in multiple human cancers and may play a role in cancer progression. Mutations in this gene are associated with the lysosomal storage disorder, Farber lipogranulomatosis, and a neuromuscular disorder, spinal muscular atrophy with progressive myoclonic epilepsy. [provided by RefSeq, Oct 2015]