

Product datasheet for RC212452

OR10A6 (NM_001004461) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OR10A6 (NM_001004461) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	OR10A6
Synonyms:	OR11-96
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC212452 representing NM_001004461 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAAGACAAAATCAAAGCTGTGTGGTTGAATTCATCCTCTTGGGCTTTTCTAACTATCCTGAGCTCC
AGGGGCAGCTCTTTGTGGCTTTCCTGGTATTTATCTGGTGACCCTGATAGGAAATGCCATTATTATAGT
CATCGTCTCCCTAGACCAGAGCCTCCACGTTCCCATGTACCTGTTTCTCCTGAACCTATCTGTGGTGGAC
CTGAGTTTCAGTGCAGTTATTATGCCTGAAATGCTGGTGGTCTCTCTACTGAAAAAACAATTTCTT
TTGGGGGCTGTTTGCACAGATGATTTTCATCCTTCTTTTGGTGGGGCTGAATGTTTTCTCTGGGAGC
AATGGCTTATGACCGATTTGCTGCAATTTGCCATCCTCTCAACTACCAAATGATTATGAATAAAGGAGTT
TTTATGAAATTAATTATATTTTCATGGGCCTTAGGTTTTATGTTAGGTACTGTTCAAACATCATGGGAT
CTAGTTTTCCCTTTGTGGCCTTAATGAAATTAACCATATATCTTGTGAAACCCAGCAGTGTTAGAACT
TGCATGTGCAGACACGTTTTTGTGTTGAAATCTATGCATTCACAGGCACCTTTTTGATTATTTGGTTCTT
TTCTTGTGTAATCTTTGCTTACATTCGAGTTCTGTTTGGCCATCCTGAAGATGCCATCAACCACTGGGA
GACAAAAGGCCTTTCCACCTGTGCCGCTCACCTCACATCTGTGACCCTATTCTATGGCAGCCAGTAT
GACTTATTTACAACCAAATCTGGCTACTACCGGAAACCAAGAAAGTGATGTATTGTCTTACTCACTT
CTGACACCACTGCTGAATCTGCTTATCTACAGTTTGCAGAAATAGTGAGATGAAGAGGGCTTTGATGAAAT
TATGGCGAAGGCGAGTGGTTTTACACACAATC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC212452 representing NM_001004461
Red=Cloning site Green=Tags(s)

MERQNQSCVVEFILLGFSNYPELQGQLFVAFVLIYLVTLIGNAIIIVIVSLDQSLHVPMYLFLNLNSVVD
 LSFSAVIMPEMLVVLSTEKTTISFGGCFAQMYFILLFPGAECFLLGAMAYDRFAAICHPLNYQMIMNKGV
 FMKLIIFSWALGFMLGTVQTSWVSSFPFCGLNEINHISCETPAVLELACADTFLFEIYAFTGTFLIILVP
 FLLILLSYIRVLFAILKMPSTTGRQKAFSTCAAHLTSVTLFYGTASMTYLQPKSGVSPETKKVMLSLSYS
 LTPLLNLLIYSLRNSEMKRALMKLWRRRVVLHTI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001004461

ORF Size: 942 bp


OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001004461.1 , NP_001004461.1
RefSeq Size:	945 bp
RefSeq ORF:	945 bp
Locus ID:	390093
UniProt ID:	Q8NH74
Cytogenetics:	11p15.4
Protein Families:	Transmembrane
Protein Pathways:	Olfactory transduction
MW:	35.1 kDa
Gene Summary:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. This olfactory receptor gene is a segregating pseudogene, where some individuals have an allele that encodes a functional olfactory receptor, while other individuals have an allele encoding a protein that is predicted to be non-functional. [provided by RefSeq, Jun 2015]

Product images:



Circular map for RC212452