

Product datasheet for **RC212469**

SPAG11B (NM_058200) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: SPAG11B (NM_058200) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: SPAG11B
Synonyms: EDDM2B; EP2; EP2C; EP2D; HE2; HE2C; SPAG11; SPAG11A
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC212469 representing NM_058200
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAGGCAACGATTGCTCCCGTCCGTACCAGCCTTCTCCTTGTTGCCCTGCTGTTCCAGGATCGTCTC
AAGCCAGACATGTGAACCACTCAGCCACTGAGGCTCTCGGAGAACTCAGGAAAGAGCCCTGGCAAGG
CACAAACGGGTTTCAGCTGCTACGCCACGCAGTAAAACGGGACCTTACCACCGCGCACCCACCTTAC
CAAGGGGATGTTCCACCGGAATTAGAAATACCATCTGCCATATGCAGCAAGGGATCTGCAGACTTTTTT
TCTGCCATTCTGGACAGGCCAGCAGCAGACAACGCTGTGGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC212469 representing NM_058200
Red=Cloning site Green=Tags(s)

MRQRLPSVTSLLLVALLPFGSSQARHVNHSATEALGELRERAPGQGTNGFQLLRHAVKRDLLPPRTPPY
QGDVPPGIRNTICHMQQGICRLLFFCHSGTGQQRHQRCC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_058200

ORF Size: 324 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_058200.2 , NP_478107.2
RefSeq Size:	755 bp
RefSeq ORF:	327 bp
Locus ID:	10407
UniProt ID:	Q08648
Cytogenetics:	8p23.1
Protein Families:	Secreted Protein
MW:	12 kDa
Gene Summary:	This gene encodes several androgen-dependent, epididymis-specific secretory proteins. The specific functions of these proteins have not been determined, but they are thought to be involved in sperm maturation. Some of the isoforms contain regions of similarity to beta-defensins, a family of antimicrobial peptides. The gene is located on chromosome 8p23 near the defensin gene cluster. Alternative splicing of this gene results in seven transcript variants encoding different isoforms. Two different N-terminal and five different C-terminal protein sequences are encoded by the splice variants. Two additional variants have been described, but their full length sequences have not been determined. [provided by RefSeq, Jul 2008]