

Product datasheet for **RC212480L1V**

SERINC3 (NM_006811) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | SERINC3 (NM_006811) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | SERINC3 |
| Synonyms: | AIGP1; DIFF33; SBBI99; TDE; TDE1; TMS-1 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_006811 |
| ORF Size: | 1419 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC212480). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_006811.2 |
| RefSeq Size: | 2592 bp |
| RefSeq ORF: | 1422 bp |
| Locus ID: | 10955 |
| UniProt ID: | Q13530 |
| Cytogenetics: | 20q13.12 |
| Domains: | TMS_TDE |
| Protein Families: | Transmembrane |



[View online »](#)

MW: 52.4 kDa

Gene Summary: Restriction factor required to restrict infectivity of lentiviruses, such as HIV-1: acts by inhibiting an early step of viral infection. Impairs the penetration of the viral particle into the cytoplasm (PubMed:26416733, PubMed:26416734).[UniProtKB/Swiss-Prot Function]