

Product datasheet for **RC212482L3V**

MS4A7 (NM_206938) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MS4A7 (NM_206938) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MS4A7
Synonyms:	4SPAN2; CD20L4; CFFM4; MS4A8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_206938
ORF Size:	585 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212482).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_206938.1 , NP_996821.1
RefSeq Size:	2889 bp
RefSeq ORF:	588 bp
Locus ID:	58475
UniProt ID:	Q9GZW8
Cytogenetics:	11q12.2
Protein Families:	Druggable Genome, Transmembrane
MW:	21.2 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the membrane-spanning 4A gene family, members of which are characterized by common structural features and similar intron/exon splice boundaries and display unique expression patterns in hematopoietic cells and nonlymphoid tissues. This family member is associated with mature cellular function in the monocytic lineage, and it may be a component of a receptor complex involved in signal transduction. This gene is localized to 11q12, in a cluster of other family members. At least four alternatively spliced transcript variants encoding two distinct isoforms have been observed. [provided by RefSeq, Jul 2008]