

### **Product datasheet for RC212501L4**

#### 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com

OriGene Technologies, Inc.

https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# ADAM12 (NM\_003474) Human Tagged Lenti ORF Clone

**Product data:** 

**Product Type:** Expression Plasmids

Product Name: ADAM12 (NM 003474) Human Tagged Lenti ORF Clone

Tag: mGFP

Symbol: ADAM12

Synonyms: ADAM12-OT1; CAR10; MCMP; MCMPMltna; MLTN; MLTNA

Mammalian Cell Puromycin

Selection:

**Vector:** pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC212501).

Sequence:

Restriction Sites: Sgfl-Mlul

**Cloning Scheme:** 





<sup>\*</sup> The last codon before the Stop codon of the ORF

**ACCN:** NM\_003474

ORF Size: 2727 bp





#### ADAM12 (NM\_003474) Human Tagged Lenti ORF Clone - RC212501L4

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

**RefSeq:** NM 003474.3, NP 003465.3

 RefSeq Size:
 5504 bp

 RefSeq ORF:
 2730 bp

 Locus ID:
 8038

 UniProt ID:
 043184

Cytogenetics:

**Domains:** Reprolysin, DISIN, Pep M12B propep, ACR

10q26.2

**Protein Families:** Druggable Genome, Protease, Secreted Protein, Transmembrane

**MW:** 99.54 kDa


**Gene Summary:** This gene encodes a member of a family of proteins that are structurally related to snake

venom disintegrins and have been implicated in a variety of biological processes involving cell-cell and cell-matrix interactions, including fertilization, muscle development, and neurogenesis. Expression of this gene has been used as a maternal serum marker for prenatal development. Alternative splicing results in multiple transcript variants encoding different isoforms. Shorter isoforms are secreted, while longer isoforms are membrane-


bound form. [provided by RefSeq, Jan 2014]



# **Product images:**



Circular map for RC212501L4



Double digestion of RC212501L4 using Sgfl and Mlul