

Product datasheet for **RC212504L4V**

SECISBP2L (NM_014701) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SECISBP2L (NM_014701) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SECISBP2L
Synonyms:	SBP2L; SLAN
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014701
ORF Size:	3168 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212504).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014701.2
RefSeq Size:	6779 bp
RefSeq ORF:	3171 bp
Locus ID:	9728
UniProt ID:	Q93073
Cytogenetics:	15q21.1
MW:	116.4 kDa



[View online »](#)

Gene Summary:

Binds SECIS (Sec insertion sequence) elements present on selenocysteine (Sec) protein mRNAs, but does not promote Sec incorporation into selenoproteins in vitro.
[UniProtKB/Swiss-Prot Function]