

## Product datasheet for **RC212520L3V**

### CD83 (NM\_001040280) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CD83 (NM_001040280) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CD83
Synonyms:	BL11; HB15
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001040280
ORF Size:	615 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212520).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001040280.1</a> , <a href="#">NP_001035370.1</a>
RefSeq Size:	2475 bp
RefSeq ORF:	615 bp
Locus ID:	9308
UniProt ID:	<a href="#">Q01151</a>
Cytogenetics:	6p23
Protein Families:	Transmembrane
MW:	23 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a single-pass type I membrane protein and member of the immunoglobulin superfamily of receptors. The encoded protein may be involved in the regulation of antigen presentation. A soluble form of this protein can bind to dendritic cells and inhibit their maturation. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011]