

Product datasheet for **RC212532L4V**

FBXO44 (NM_001014765) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FBXO44 (NM_001014765) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FBXO44
Synonyms:	FBG3; FBX6A; FBX30; Fbx44; Fbxo6a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001014765
ORF Size:	765 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212532).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001014765.1 , NP_001014765.1
RefSeq Size:	3320 bp
RefSeq ORF:	768 bp
Locus ID:	93611
UniProt ID:	Q9H4M3
Cytogenetics:	1p36.22
Protein Families:	Druggable Genome
MW:	29.6 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the F-box protein family which is characterized by an approximately 40 amino acid motif, the F-box. The F-box proteins constitute one of the four subunits of the ubiquitin protein ligase complex called SCFs (SKP1-cullin-F-box), which function in phosphorylation-dependent ubiquitination. The F-box proteins are divided into 3 classes: Fbws containing WD-40 domains, Fbls containing leucine-rich repeats, and Fbxs containing either different protein-protein interaction modules or no recognizable motifs. The protein encoded by this gene belongs to the Fbxs class. It is also a member of the NFB42 (neural F Box 42 kDa) family, similar to F-box only protein 2 and F-box only protein 6. Several alternatively spliced transcript variants encoding two distinct isoforms have been found for this gene. [provided by RefSeq, Feb 2015]