

## Product datasheet for **RC212535L4V**

### Zfp219 (ZNF219) (NM\_001102454) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Zfp219 (ZNF219) (NM_001102454) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Zfp219
Synonyms:	ZFP219
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001102454
ORF Size:	1389 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212535).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001102454.1</a> , <a href="#">NP_001095924.1</a>
RefSeq Size:	2954 bp
RefSeq ORF:	2169 bp
Locus ID:	51222
UniProt ID:	<a href="#">Q9P2Y4</a>
Cytogenetics:	14q11.2
Protein Families:	Transcription Factors
MW:	76.7 kDa



[View online »](#)

**Gene Summary:**

This gene is a member of the Kruppel-like zinc finger gene family. The encoded protein functions as a transcriptional repressor of the high mobility group nucleosome binding domain 1 protein, which is associated with transcriptionally active chromatin. [provided by RefSeq, Apr 2017]