

Product datasheet for **RC212563A1V**

Human TNNT3 (NM_001042780) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TNNT3 (NM_001042780) AAV Particle
Tag:	Myc-DDK
Symbol:	TNNT3
Synonyms:	beta-TnTF; DA2B2; TNTF
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC212563 representing NM_001042780 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGCTGACGAGGAAGTTGAACAGGTGGAGGAGCAGTACGAAGAAGAAGAGGAAGCCAGGAGGAAGAGG
AAGTTCAAGAAGAGGAGAAACCGAGACCCAACTCACTGCTCCTAAGATCCCAGAAGGGGAGAAAGTGGA
CTTCGATGACATCCAGAAGAAGCGTCAGAACAAAGACCTAATGGAGCTCCAGGCCCTCATCGACAGCCAC
TTTGAAGCCCGGAAGAAGGAGGAGGAGGAGCTGGTCGCTCTCAAAGAGAGAATCGAGAAGCGCCGTCGAG
AGAGAGCGGAGCAGCAGAGGATTTCGTGCAGAGAAGGAGAGGGAGCGCCAGAACAGACTGGCGGAGGAAAA
GGCCAGAAGGGAGGAGGAGGATGCCAAGAGGAGGGCAGAGGACGACCTGAAGAAGAAGAAAGCTCTGTCT
TCCATGGGAGCCAACTACAGCAGCTACCTGGCCAAGGCTGACCAGAAGAGAGGCAAGAAGCAGACAGCCC
GGGAAATGAAGAAGAAGATTCTGGCTGAGAGACGCAAGCCGCTCAACATCGATCACCTTGGTGAAGACAA
ACTGAGGGACAAGGCCAAGGAGCTCTGGGAGACCCTGCACCAGCTGGAGATTGACAAGTTCGAGTTTGGG
GAGAAGCTGAAACGCCAGAAATATGACATCATGAATGTCCGGGCCAGAGTGCAGATGCTGGCCAAGTTCA
GCAAGAAGGCTGGGACCCAGCCAAGGGCAAAGTCGGCGGGCGCTGGAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC212563 representing NM_001042780 Red=Cloning site Green=Tags(s) MSDEEVEQVEEQYEEEEEAQEEEEVQEEKPRPKLTAPKIPGEKVDFFDIQKKRQNKDLMELQALIDSH FEARKKEEEELVALKERIEKRRRAEAEQQRIRAEKERERQNLAEKARREEDAKRRAEDDLKKKKALS SMGANYSSYLAKADQKRGKKQTAREMKKKILAERRKPLNIDHLGEDKLRDKAKELWETLHQLEIDKFEFG EKLKRQKYDIMNVRARVQMLAKFSKKAGTPAKGKVGGRWK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001042780
ORF Size:	750 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001042780.1</u>
RefSeq Size:	1193 bp
RefSeq ORF:	753 bp
Locus ID:	7140
UniProt ID:	<u>P45378</u>
Cytogenetics:	11p15.5
MW:	30.2 kDa