

Product datasheet for RC212573

OR4F4 (NM_001004195) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	OR4F4 (NM_001004195) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	OR4F4
Synonyms:	OLA-7501; OR4F18
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC212573 representing NM_001004195 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTGACTGAATTCATTTTTCTGGGTCTCTGATTCTCAGGAACCTCCAGACCTTCTATTTATGTTGT
TTTTTGATTCTATGGAGGAATCGTGTTGGAAACCTTCTATTGTCATAACAGTGGTATCTGACTCCCA
CCTTCACTCTCCCATGTACTTCTGCTAGCCAACCTCTCACTCATTGATCTGTCTGTCTTCAGTCACA
GCCCCAAGATGATTACTGACTTTTTAGCCAGCGCAAAGTCATCTTTCAAGGGCTGCCTTGTTCAGA
TATTTCTCCTTCACTTCTTTGGTGGGAGTGAGATGGTGATCCCTCATAGCCATGGGCTTTGACAGATAT
AGCAATATGCAAACCCCTACACTACACTACAATTATGTGTGGCAACGCATGTGTCCGCATTATGGCTGTC
GCATGGGGAATTGGCTTTCTCCATTCCGGTGAGCCAGTTGGCCTTTGCCGTGCACTTACCCTTCTGTGGTC
CCAATGAGGTCGATAGTTTTTATTGTGACCTTCTAGGGTAATCAAACCTGCCTGTACAGATACCTACAG
GCTAGATATTATGGTCATTGCTAACAGTGGTGTGCTCACTGTGTGTTCTTTGTTCTTCTAATCATCTCA
TACACTATCATCCTAATGACCATCCAGCATTGCCCTTAGATAAGTCGTCCAAAGCTCTGTCCACTTTGA
CTGCTCACATTACAGTAGTTCTTTGTTCTTTGGACCATGTGTCTTTATTTATGCCTGGCCATCCCCAT
CAAGTCATTAGATAAATTCCTTGCTGTATTTATTCTGTGATCACCCCTCTTTGAACCAATTATATAC
ACACTGAGGAACAAAGACATGAAGACGGCAATAAGACGGCTGAGAAAATGGGATGCACATTCTAGTGTA
AGTTT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC212573 representing NM_001004195
Red=Cloning site Green=Tags(s)

MVTEFIFLGLSDSQELQTFLFMLFFVFYGGIVFGNLLIVITVVSDSHLHSPMYFLLANLSLIDLSLSSVT
 APKMITDFFSQRKVISFKGCLVQIFLLHFFGGSEMVLIAMGFDRYIAICKPLHYTTIMCGNACVGIMAV
 AWGIGFLHSVSQLAFVHLPCGPNVDSFYCDLPRVIKLACTDTYRLDIMVIANSGLVTVCSFVLLIIS
 YTIILMTIQHCPLDKSSKALSTLTAHITVVLLFFGPCVFIYAWPFPIKSLDKFLAVFYSVITPLLNPDIY
 TLRNKDMKTAIRRLRKWDAHSSVKF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8001_h04.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001004195

ORF Size: 915 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001004195.2](#), [NP_001004195.2](#)

RefSeq Size: 918 bp

RefSeq ORF: 918 bp

Locus ID: 26682

UniProt ID: [Q96R69](#)

Cytogenetics: 15q26.3


Protein Families: Transmembrane

Protein Pathways: Olfactory transduction

MW: 34.1 kDa

Gene Summary: Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC212573