

Product datasheet for **RC212583L3V**

Her2 (ERBB2) (NM_004448) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Her2 (ERBB2) (NM_004448) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Her2
Synonyms:	CD340; HER-2; HER-2/neu; HER2; MLN 19; NEU; NGL; TKR1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004448
ORF Size:	3765 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212583).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004448.2
RefSeq Size:	4624 bp
RefSeq ORF:	3768 bp
Locus ID:	2064
UniProt ID:	P04626
Cytogenetics:	17q12
Domains:	Recep_L_domain, pkinase, TyrKc, S_TKc, YLP, Furin-like, FU
Protein Families:	Druggable Genome, Protein Kinase, Transmembrane




[View online »](#)

Protein Pathways: Adherens junction, Bladder cancer, Calcium signaling pathway, Endometrial cancer, ErbB signaling pathway, Focal adhesion, Non-small cell lung cancer, Pancreatic cancer, Pathways in cancer, Prostate cancer


MW: 137.7 kDa

Gene Summary: This gene encodes a member of the epidermal growth factor (EGF) receptor family of receptor tyrosine kinases. This protein has no ligand binding domain of its own and therefore cannot bind growth factors. However, it does bind tightly to other ligand-bound EGF receptor family members to form a heterodimer, stabilizing ligand binding and enhancing kinase-mediated activation of downstream signalling pathways, such as those involving mitogen-activated protein kinase and phosphatidylinositol-3 kinase. Allelic variations at amino acid positions 654 and 655 of isoform a (positions 624 and 625 of isoform b) have been reported, with the most common allele, Ile654/Ile655, shown here. Amplification and/or overexpression of this gene has been reported in numerous cancers, including breast and ovarian tumors. Alternative splicing results in several additional transcript variants, some encoding different isoforms and others that have not been fully characterized. [provided by RefSeq, Jul 2008]

Product images:



[RC212583L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC212583L3V particle to overexpress human ERBB2-Myc-DDK fusion protein.



[RC212583L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC212583L3V particle to overexpress human ERBB2-Myc-DDK fusion protein.