

Product datasheet for **RC212586L2V**

TPPP3 (NM_015964) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TPPP3 (NM_015964) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TPPP3
Synonyms:	CGI-38; p20; p25gamma; TPPP/p20
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_015964
ORF Size:	528 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212586).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015964.2
RefSeq Size:	1055 bp
RefSeq ORF:	531 bp
Locus ID:	51673
UniProt ID:	Q9BW30
Cytogenetics:	16q22.1
MW:	18.8 kDa


[View online »](#)

Gene Summary: Binds tubulin and has microtubule bundling activity. May play a role in cell proliferation and mitosis.[UniProtKB/Swiss-Prot Function]