

Product datasheet for **RC212596A1V**

Human SNAP25 (NM_130811) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human SNAP25 (NM_130811) AAV Particle
Tag:	Myc-DDK
Symbol:	SNAP25
Synonyms:	bA416N4.2; CMS18; dj1068F16.2; RIC-4; RIC4; SEC9; SNAP; SNAP-25; SUP
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC212596 representing NM_130811 Red=Cloning site Blue=ORF Green=Tags(s) TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC ATGGCCGAAGACGCAGACATGCGCAATGAGCTGGAGGAGATGCAGCGAAGGGCTGACCAGTTGGCTGATGAGTCGCTGGAAGCACCCTCGTATGCTGCAACTGGTTGAAGAGAGTAAAGATGCTGGTATCAGGACTTTGGTTATGTTGGATGAACAAGGAGAACAACCTGGAACGCATTGAGGAAGGGATGGACCAATCAATAAGGACATGAAAGAAGCAGAAAAGAATTTGACGGACCTAGGAAAATTCGCGGGCTTTGTGTGTCCCTGTAACAAGCTTAAATCAAGTGATGCTTACAAAAAGCCTGGGGCAATAATCAGGACGGAGTGGTGGCCAGCCAGCCTGCTCGTGTAGTGGACGAACGGGAGCAGATGGCCATCAGTGGCGGCTTCATCCGCAGGGTAACAAATGATGCCCGAGAAAATGAAATGGATGAAAACCTAGAGCAGGTGAGCGGCATCATCGGGAACCTCCGTACATGGCCCTGGATATGGGCAATGAGATCGATACACAGAATCGCCAGATCGACAGGATCATGGAGAAGGCTGATTCACACAAAACCAGAATTGATGAGGCCAACCAACGTGCAACAAAGATGCTGGGAAGTGGT ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
Protein Sequence:	>RC212596 representing NM_130811 Red=Cloning site Green=Tags(s) MAEDADMRNELEEMQRRADQLADESLESTRMLQLVEESKDAGIRTLVMLDEQGEQLERIEEGMDQINKDMKEAEKNLTDLGKFCGLCVCPCNKLKSSDAYKKAWGNNQDGVVASQPARVVDEREQMAISGGFIRRVNTDARENEMDENLEQVSGIIGNLRHMLDMGNEIDTQNRQIDRIMEKADSNKTRIDEANQRATKMLGSG TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_130811
ORF Size:	618 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_130811.1</u>
RefSeq Size:	2053 bp
RefSeq ORF:	621 bp
Locus ID:	6616
UniProt ID:	<u>P60880</u>
Cytogenetics:	20p12.2
MW:	23.1 kDa