

## Product datasheet for **RC212596L3V**

### SNAP25 (NM\_130811) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | SNAP25 (NM_130811) Human Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | SNAP25   |
| Synonyms:                 | bA416N4.2; CMS18; dj1068F16.2; RIC-4; RIC4; SEC9; SNAP; SNAP-25; SUP   |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)   |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_130811  |
| ORF Size:                 | 618 bp   |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC212596).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_130811.1</a>  |
| RefSeq Size:              | 2053 bp  |
| RefSeq ORF:               | 621 bp   |
| Locus ID:                 | 6616   |
| UniProt ID:               | <a href="#">P60880</a>   |
| Cytogenetics:             | 20p12.2  |
| Domains:                  | t_SNARE, SNAP-25   |
| Protein Families:         | Druggable Genome   |



[View online »](#)

**Protein Pathways:** SNARE interactions in vesicular transport

**MW:** 23.1 kDa

**Gene Summary:** Synaptic vesicle membrane docking and fusion is mediated by SNAREs (soluble N-ethylmaleimide-sensitive factor attachment protein receptors) located on the vesicle membrane (v-SNAREs) and the target membrane (t-SNAREs). The assembled v-SNARE/t-SNARE complex consists of a bundle of four helices, one of which is supplied by v-SNARE and the other three by t-SNARE. For t-SNAREs on the plasma membrane, the protein syntaxin supplies one helix and the protein encoded by this gene contributes the other two. Therefore, this gene product is a presynaptic plasma membrane protein involved in the regulation of neurotransmitter release. Two alternative transcript variants encoding different protein isoforms have been described for this gene. [provided by RefSeq, Jul 2008]