

Product datasheet for **RC212604L3V**

LLGL2 (NM_004524) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	LLGL2 (NM_004524) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LLGL2
Synonyms:	HGL; Hugl-2; LGL2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004524
ORF Size:	3045 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212604).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004524.2
RefSeq Size:	3547 bp
RefSeq ORF:	3048 bp
Locus ID:	3993
UniProt ID:	Q6P1M3
Cytogenetics:	17q25.1
Domains:	WD40
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Tight junction

MW: 113.4 kDa

Gene Summary: The lethal (2) giant larvae protein of *Drosophila* plays a role in asymmetric cell division, epithelial cell polarity, and cell migration. This human gene encodes a protein similar to lethal (2) giant larvae of *Drosophila*. In fly, the protein's ability to localize cell fate determinants is regulated by the atypical protein kinase C (aPKC). In human, this protein interacts with aPKC-containing complexes and is cortically localized in mitotic cells. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]