

## Product datasheet for **RC212615L3V**

### **PAX3 (NM\_013942) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	PAX3 (NM_013942) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PAX3
Synonyms:	CDHS; HUP2; WS1; WS3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_013942
ORF Size:	618 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212615).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_013942.4</a>
RefSeq Size:	1095 bp
RefSeq ORF:	621 bp
Locus ID:	5077
UniProt ID:	<a href="#">P23760</a>
Cytogenetics:	2q36.1
Protein Families:	Adult stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Transcription Factors




[View online »](#)

**MW:** 22.7 kDa

**Gene Summary:** This gene is a member of the paired box (PAX) family of transcription factors. Members of the PAX family typically contain a paired box domain and a paired-type homeodomain. These genes play critical roles during fetal development. Mutations in paired box gene 3 are associated with Waardenburg syndrome, craniofacial-deafness-hand syndrome, and alveolar rhabdomyosarcoma. The translocation t(2;13)(q35;q14), which represents a fusion between PAX3 and the forkhead gene, is a frequent finding in alveolar rhabdomyosarcoma. Alternative splicing results in transcripts encoding isoforms with different C-termini. [provided by RefSeq, Jul 2008]

### Product images:



[RC212615L3] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC212615L3V particle to overexpress human PAX3-Myc-DDK fusion protein.