

## Product datasheet for **RC212639L3V**

### SNRK (NM\_001100594) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SNRK (NM_001100594) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SNRK
Synonyms:	HSNFRK
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001100594
ORF Size:	2295 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212639).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001100594.1</a> , <a href="#">NP_001094064.1</a>
RefSeq Size:	5174 bp
RefSeq ORF:	2298 bp
Locus ID:	54861
UniProt ID:	<a href="#">Q9NRH2</a>
Cytogenetics:	3p22.1
Protein Families:	Druggable Genome, Protein Kinase
MW:	84.4 kDa



[View online »](#)

**Gene Summary:**

SNRK is a member of the sucrose nonfermenting (SNF)-related kinase family of serine/threonine kinases (Kertesz et al., 2002 [PubMed 12234663]).[supplied by OMIM, Apr 2009]