

Product datasheet for **RC212656A1V**

Human PDLIM7 (NM_203352) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human PDLIM7 (NM_203352) AAV Particle
Tag:	Myc-DDK
Symbol:	PDLIM7
Synonyms:	LMP1; LMP3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC212656 representing NM_203352 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGATTCTTCAAAGTAGTGCTGGAGGGGCCAGCACCTTGGGGCTTCCGGCTGCAAGGGGGCAAGGACT
TCAATGTGCCCCCTCTCCATTTCCCGGCTCACTCTGGGGGCAAAGCGGCGCAGGCCGGAGTGCCCGTGGG
TGACTGGGTGCTGAGCATCGATGGCGAGAATGCGGGTAGCCTCACACACATCGAAGCTCAGAACAAAGATC
CGGGCCTGCGGGGAGCGCCTCAGCCTGGGCCTCAGCAGGGCCAGCCGGTTCAGAGCAAACCGCAGAAGG
TGCAGACCCCTGACAAACAGCCGCTCCGACCGCTGGTCCCAGATGCCAGCAAGCAGCGGCTGATGGAGAA
CACAGAGGACTGGCGGGCCGGCCGGGACAGGCCAGTCGCGTTCCCTCCGCATCCTTGCCACCTCACAC
GGCACCGAGTTTCATGCAAGACCCGGATGAGGAGCACCTGAAGAAATCAAGCCAGGTGCCAGGACAGAAG
CCCCAGCCCCAGCCTCATCTACACCCAGGAGCCCTGGCCTGGCCCTACCGCCCCAGCCCTACCAGCCG
CCCGCCTGGGCTGTGGACCCTGCGTTTGCCGAGCGCTATGCCCGGACAAAACGAGCACAGTGCTGACC
CGGCACAGCCAGCCGGCCACGCCACGCCGCTGCAGAGCCGCACCTCCATTGTGCAGGCAGCTGCCGGAG
GGGTGCCAGGAGGGGGCAGCAACAACGGCAAGACTCCCGTGTGTACCAGTGCCACAAGGTATCCGGGG
CCGCTACCTGGTGGCGCTGGGCCACGCGTACCACCGGAGGAGTTTGTGTGTAGCCAGTGTGGGAAGGTC
CTGGAAGAGGGTGGCTTCTTTGAGGAGAAGGGGCCATCTTCTGCCACCATGCTATGACGTGCGCTATG
CACCCAGCTGTGCAAGTGCAAGAAGAAGATTACAGGCGAGATCATGCACGCCCTGAAGATGACCTGGCA
CGTGCACTGCTTTACCTGTGCTGCCTGCAAGACGCCCATCCGGAACAGGGCCTTCTACATGGAGGAGGGC
GTGCCCTATTGCGAGCGAGACTATGAGAAGATGTTTGGCACGAAATGCCATGGCTGTGACTTCAAGATCG
ACGCTGGGGACCGCTTCTGGAGGCCCTGGGCTTCAGCTGGCATGACACCTGCTTCGTCTGTGCGATATG
TCAGATCAACCTGGAAGGAAAGACCTTCTACTCCAAGAAGGACAGGCCTCTCTGCAAGAGCCATGCCTTC
TCTCATGTG

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC212656 representing NM_203352 Red=Cloning site Green=Tags(s) MDSFKVVLGPAPWGFRLQGGKDFNVPLSISRLTPGGKAAQAGVAVGDWVLSIDGENAGSLTHIEAQNKI RACGERLSLGLSRAQPVQSKPQKVQTPDKQPLRPLVPDASKQRLMENTEDWRPRPGTGQSRSFRLAHLT GTEFMQDPDEEHLKSSQVPRTEAPAPASSTPQEPWPGPTAPSPTRPPWAVDPFAAERYAPDKTSTVL RHSQPATPTPLQSRSTIVQAAAGGVPGGGSNNGKTPVCHQCHKVIRGRYLVALGHAYHPEEFVCSQCGKV LEEGGFEEKGAIFCPPCYDVRYAPSCAKCKKITGEIMHALKMTWHVHCFTCAACKTPIRNRAFYMEEG VPYCERDYEKMFGTKCHGCDFKIDAGDRFLEALGFSWHDTCFVCAICQINLEGKTFYSKKDRPLCKSHAF SHV SGPTRTRRLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_203352
ORF Size:	1269 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_203352.1
RefSeq Size:	1607 bp
RefSeq ORF:	1272 bp
Locus ID:	9260
UniProt ID:	Q9NR12
Cytogenetics:	5q35.3
MW:	46.3 kDa