

## Product datasheet for **RC212672L4V**

### **RASGRP2 (NM\_001098670) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	RASGRP2 (NM_001098670) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RASGRP2
Synonyms:	CALDAG-GEFI; CDC25L
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001098670
ORF Size:	1827 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212672).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001098670.1</a>
RefSeq Size:	2244 bp
RefSeq ORF:	1830 bp
Locus ID:	10235
UniProt ID:	<a href="#">Q7LDG7</a>
Cytogenetics:	11q13.1
Protein Pathways:	Chemokine signaling pathway, MAPK signaling pathway
MW:	69.2 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a brain-enriched nucleotide exchanged factor that contains an N-terminal GEF domain, 2 tandem repeats of EF-hand calcium-binding motifs, and a C-terminal diacylglycerol/phorbol ester-binding domain. This protein can activate small GTPases, including RAS and RAP1/RAS3. The nucleotide exchange activity of this protein can be stimulated by calcium and diacylglycerol. Four alternatively spliced transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Jan 2016]