

Product datasheet for **RC212692L2V**

NFAT2 (NFATC1) (NM_006162) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | NFAT2 (NFATC1) (NM_006162) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | NFATC1 |
| Synonyms: | NF-ATC; NF-ATc1.2; NFAT2; NFATc |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-mGFP (PS100071) |
| Tag: | mGFP |
| ACCN: | NM_006162 |
| ORF Size: | 2475 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC212692). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_006162.3 |
| RefSeq Size: | 4762 bp |
| RefSeq ORF: | 2478 bp |
| Locus ID: | 4772 |
| UniProt ID: | O95644 |
| Cytogenetics: | 18q23 |
| Domains: | RHD, IPT |
| Protein Families: | Druggable Genome, Transcription Factors |



[View online »](#)

| | |
|--------------------------|---|
| Protein Pathways: | Axon guidance, B cell receptor signaling pathway, Natural killer cell mediated cytotoxicity, T cell receptor signaling pathway, VEGF signaling pathway, Wnt signaling pathway |
| MW: | 88.6 kDa |
| Gene Summary: | <p>The product of this gene is a component of the nuclear factor of activated T cells DNA-binding transcription complex. This complex consists of at least two components: a preexisting cytosolic component that translocates to the nucleus upon T cell receptor (TCR) stimulation, and an inducible nuclear component. Proteins belonging to this family of transcription factors play a central role in inducible gene transcription during immune response. The product of this gene is an inducible nuclear component. It functions as a major molecular target for the immunosuppressive drugs such as cyclosporin A. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene. Different isoforms of this protein may regulate inducible expression of different cytokine genes. [provided by RefSeq, Jul 2013]</p> |