

Product datasheet for **RC212737A1V**

Human DcR3 (TNFRSF6B) (NM_032945) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human DcR3 (TNFRSF6B) (NM_032945) AAV Particle
Tag:	Myc-DDK
Symbol:	DcR3
Synonyms:	DCR3; decoy receptor 3; DJ583P15.1.1; M68; OTTHUMP00000031583; TR6; tumor necrosis factor receptor superfamily, member 6b; tumor necrosis factor receptor superfamily, member 6b, decoy
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC212737 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGGGCGCTGGAGGGGCCAGGCCTGTCGCTGCTGTGCCTGGTGTGGCGCTGCCTGCCCTGCTGCCGG
TGCCGGCTGTACGCGGAGTGGCAGAAACACCCACCTACCCCTGGCGGGACGCAGAGACAGGGGAGCGGCT
GGTGTGCGCCCCAGTGCCCCCAGGCACCTTTGTGCAGCGGCCGTGCCGCCGAGACAGCCCCACGACGTGT
GGCCCGTGTCCACCGCGCCACTACACGCAGTTCTGGAACCTACCTGGAGCGCTGCCGCTACTGCAACGTCC
TCTGCGGGGAGCGTGAGGAGGAGGCACGGGCTTGCCACGCCACCCACAACCGTGCCTGCCGCTGCCGCAC
CGGCTTCTTCGCGCACGCTGGTTTCTGCTTGAGCACGCATCGTGTCCACCTGGTGCCGGCGTGATTGCC
CCGGGCACCCCGCCAGAACACGCAGTGCCAGCCGTGCCCCCAGGCACCTTCTCAGCCAGCAGCTCCA
GCTCAGAGCAGTGCCAGCCCCACGCAACTGCACGGCCCTGGGCTGGCCCTCAATGTGCCAGGCTCTTC
CTCCCATGACACCCTGTGCACCAGCTGCACTGGCTTCCCCCTCAGCACCAGGTTACCAGGAGCTGAGGAG
TGTGAGCGTGCCGTCATCGACTTTGTGGCTTTCCAGGACATCTCCATCAAGAGGCTGCAGCGGCTGCTGC
AGGCCCTCGAGGCCCGGAGGGCTGGGTCCGACACCAAGGGCGGGCCGCGCGGCTTGCAGCTGAAGCT
GCGTCGGCGGCTCACGGAGCTCCTGGGGCGCAGGACGGGGCGCTGCTGGTGCGGCTGCTGCAGGCGCTG
CGCGTGCCAGGATGCCCGGGCTGGAGCGGAGCGTCCGTGAGCGCTTCTCCCTGTGCAC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC212737 protein sequence Red=Cloning site Green=Tags(s) MRALEGPGLSLLCLVLALPALLPVPVAVRGVAETPTYPWRDAETGERLVCAQCPPGTFVQRPCRRDSPTTC GPCPPRHYTQFWNYLERCRYCNVLCGEREEARACHATHNRACRCRTGFFAHAGFCLEHASCPPGAGVIA PGTPSQNTQCQPCPPGTFSSSSSEQCQPHRNCTALGLALNVPGSSSHDTLCTSCTGFPLSTRVPGAEE CERAVIDFVAFQDISIKRLQRLQALEAPEGWGPTPRAGRAALQLKLRRRLTELLGAQDGALLVRLQAL RVARMPGLERSVRERFLPVH TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_032945
ORF Size:	900 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_032945.2 , NP_116563.1
RefSeq Size:	1458 bp
RefSeq ORF:	902 bp
Locus ID:	8771
Cytogenetics:	20q13.33
MW:	32.7 kDa