

Product datasheet for **RC212749**

Invadolysin (LMLN) (NM_033029) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Invadolysin (LMLN) (NM_033029) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	LMLN
Synonyms:	GP63; INV; IX14; LMNL1; MSP
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RC212749 representing NM_033029
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGTAACGACGCTCGGCCGAAGATGGCGGCCGAATGGGGCGGAGGAGTGGGTTACTCGGGCTCAGGCC
 CGGGCCGAGCCGGTGGCGCTGGAGCGGGTCTGTGTGGTCCGAAGCGTTTTACTCCTGTTGGCGGGCT
 CCGGGCCAGCGCCACATCTACTCCCGTCTCCTTGGGCAGTTCCTCCCTGCCCGCACCAGTCCCTCT
 GACACTGAGGTCATAAATAAAGTTCATCTTAAGGCAAATCATGTGGTCAAGAGAGATGTTGATGAGCATT
 TAAGAATCAAGACTGTCTATGATAAAAGTGTGAAGAGTTGCTCCCTGAGAAAAAGAATCTGTAAAGAA
 CAAGCTTTTCCACAAGCGATTTCTTATTTAGAGAAGACTTTTCAGGTCCGTCGACTGCGGGCACTATC
 TTACTIONTAGCAGACAATGTGCAACAAACCAATACCTCCGGAAGGAAAACGATCCTCACAGGTAAGTGCACCG
 GGGAGTGTGCCGACACACAAAGTGGCGCCCGTTATTGTTCTGAGGAACATCTCCAGCAATGCCGGT
 CTACCGTGGGGTAAGTGGCCTCATGGAGCAGTGGGTGTCCAGACCAAGAAGGCATCTCAGATGCAGAC
 TTTGTTCTTTACGTTGGTCTCTGGCCACCGAGAGATGCAGCCATGAAAACATCATCTCTTATGCAGCCT
 ATTGTCAGCAGGAAGCAAACATGGACAGGCCAATAGCAGGATATGCTAACCTGTGTCCAAATATGATCTC
 TACCCAGCCTCAGGAGTTTGTGGGATGCTGTCCACAGTAAAACATGAGGTTATTCATGCCTGGGTTTC
 TCTGCTGGGCTGTTTGCATTCTACCATGATAAAGATGAAATCCTCTCACTTCAAGATTTGCAGATGGCC
 TCCCACCTTTAATTATAGTCTGGGATTATCAATGGAGTGATAAAGTAGTTCGAAAAGTGGAGAGATT
 ATGGGATGTTTCGAGATAAAGATAGTTCGTCACACTGTGTATCTCCTGGTAACGCCCTCGTGTGTTGAG
 GAAGCACGAAAACATTTGATTGTCCAGTCTAGAGGGAATGAACTTAAAAATCAAGGTGGTGTGGGCA
 CTGAGCTCAACCATTTGGGAAAAAAGGTTATTAGAGAATGAAGCGATGACTGGTTCTCACACTCAGAATCG
 AGTACTCTCTCGAATCACTCTGGCATTAAATGGAGGACACTGGGAGACAGATGCTGAGCCCTTACTGTGAC
 ACGCTCAGAAGTAACCCACTGCAGCTAAGTTCAGACAGGACCAGAGAGCAGTTGCCGTGTGAATTTGC
 AGAAGTCCCTAAGCCTTACCACAGGAATACCAGTACTTTGATGAACTCAGTGAATACCTGCAGAAGA
 TTTGCCTTATTATGGTGGCTCCGTGAAATGCTGACTACTGCCCTTTCAGTCAGGAATTCAGTTGGCAT
 TTAAGTGGTGAATATCAGCGCAGCTCAGATTGTAGAATATTGGAAAATCAACCAGAAAATTTTAAGAACT
 ATGGCGCTGAAAAGTATGGACCTCATTCCGTTTGTCTAATTCAGAAATCAGCATTCTGTTATGGAGAAGTG
 TGAGAGGAAGCTGAGTTACCCAGACTGGGAAGCGGATGCTATCAGGTTTCTTGTCTCCTCAAGGTCTG
 AAAGTTTGGGTCCAAGATACTTCATATTTGTGTAGTCGGGCTGGGCAGGTCTCCCTGTGAGTATCCAGA
 TGAAATGGCTGGATTCACGATGAAAACCTGCTCTGCCCATCATGTTGGGACTTCTGTGAGCTGTCTCC
 AGAAACAGATCCTCCAGCCACTAACCTGACCCGAGCTCTGCCACTTGATCTTTGTTCTGTTCTCGAGC
 CTGGTGGTCACCCCTCTGGCTTCTGCTAGGCAATCTGTTTCTCTGCTGGCTGGATTTCTTCTGTGTATAT
 GGCAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC212749 representing NM_033029
 Red=Cloning site Green=Tags(s)

```

MVTTLGPKMAAEWGGGVGYSGSGPGRSRWRWVRSVLLLLGGLRASATSTPVSLGSSPPCRHHVPS
DTEVINKVHLKANHVVKRDVDEHLRIKTVYDKSVEELLPEKKNLVKNLFPQAISYLEKTFQVRRPAGTI
LLSRQCATNQYLKENDPHRYCTGECAAHTKCGPVIVPEEHLQQCRVYRGGKWPBGAVGVPDQEGISDAD
FVLYVYGALATERCSHENIISYAAYCQEQEANMDRPIAGYANLCPNMISTQPQEFVGMLSTVKHEVIHALGF
SAGLFAFYHDKDGNPLTSRFADGLPPFNYSGLYQWSDKVVRKVERLWDVRDNKIVRHTVYLLVTPRVVE
EARKHFDCPVLEGMELNQGGVGTLELNHWKRLLENEAMTGSHTQNRVLSRITLALMEDTGRQMLSPYCD
TLRSNPLQLTCRQDQRAVAVCNLQKFKPLPQEQYFDELSPGPAEDLPYYGGSVEIADYCPFSQEFSWH
LSGEYQRSSDCRILENQPEIFKNYGAEKYGPSVCLIQKSAFVMEK CERKLSYPDWGSGCYQVSCSPQGL
KVWVQDTSYLC SRAGQVLPVSIQMNGWIHDGNLLCPSCWDFCELCPPETDPPATNLTRALPLDLCSCSSS
LVVTLWLLLGNLFPLLAGFLLCIWH
  
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_033029

ORF Size: 1965 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_033029.3](#), [NP_149018.2](#)

RefSeq Size: 7043 bp

RefSeq ORF: 1944 bp

Locus ID: 89782

UniProt ID: [Q96KR4](#)

Cytogenetics: 3q29

Protein Families: Druggable Genome, Protease, Transmembrane

MW: 73.6 kDa

Gene Summary: This gene encodes a zinc-metallopeptidase. The encoded protein may play a role in cell migration and invasion. Studies of a similar protein in *Drosophila* indicate a potential role in mitotic progression. Alternatively spliced transcript variants have been described. [provided by RefSeq, Feb 2009]