

Product datasheet for **RC212750L3V**

NKAIN2 (NM_001040214) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NKAIN2 (NM_001040214) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NKAIN2
Synonyms:	FAM77B; NKAIP2; TCBA; TCBA1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001040214
ORF Size:	624 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212750).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001040214.1
RefSeq Size:	3315 bp
RefSeq ORF:	627 bp
Locus ID:	154215
UniProt ID:	Q5VXU1
Cytogenetics:	6q22.31
Protein Families:	Transmembrane
MW:	23.7 kDa



[View online »](#)

Gene Summary:

This gene encodes a transmembrane protein that interacts with the beta subunit of a sodium/potassium-transporting ATPase. A chromosomal translocation involving this gene is a cause of lymphoma. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2014]