

## Product datasheet for **RC212769L3V**

### D aspartate oxidase (DDO) (NM\_004032) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	D aspartate oxidase (DDO) (NM_004032) Human Tagged ORF Clone Lentiviral Particle
Symbol:	D aspartate oxidase
Synonyms:	DASOX; DDO-1; DDO-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004032
ORF Size:	930 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212769).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004032.2</a> , <a href="#">NP_004023.2</a>
RefSeq Size:	1556 bp
RefSeq ORF:	849 bp
Locus ID:	8528
UniProt ID:	<a href="#">Q99489</a>
Cytogenetics:	6q21
Domains:	DAO
Protein Pathways:	Alanine, aspartate and glutamate metabolism



[View online »](#)

**MW:** 33.8 kDa

**Gene Summary:** The protein encoded by this gene is a peroxisomal flavoprotein that catalyzes the oxidative deamination of D-aspartate and N-methyl D-aspartate. Flavin adenine dinucleotide or 6-hydroxyflavin adenine dinucleotide can serve as the cofactor in this reaction. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2019]