

Product datasheet for **RC212770L4V**

USP32 (NM_032582) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	USP32 (NM_032582) Human Tagged ORF Clone Lentiviral Particle
Symbol:	USP32
Synonyms:	NY-REN-60; USP10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_032582
ORF Size:	4811 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212770).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_032582.3
RefSeq Size:	7026 bp
RefSeq ORF:	4815 bp
Locus ID:	84669
UniProt ID:	Q8NFA0
Cytogenetics:	17q23.1-q23.2
Domains:	UCH, EFh, DUSP
Protein Families:	Druggable Genome, Protease



[View online »](#)

MW: 181.5 kDa