

## Product datasheet for **RC212806L4V**

### THSD1 (NM\_018676) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	THSD1 (NM_018676) Human Tagged ORF Clone Lentiviral Particle
Symbol:	THSD1
Synonyms:	ANIB12; TMTSP; UNQ3010
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018676
ORF Size:	2556 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212806).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_018676.2</a>
RefSeq Size:	3054 bp
RefSeq ORF:	2559 bp
Locus ID:	55901
UniProt ID:	<a href="#">Q9NS62</a>
Cytogenetics:	13q14.3
Domains:	tsp_1
Protein Families:	Transmembrane


[View online »](#)

**MW:** 94.4 kDa

**Gene Summary:** The protein encoded by this gene contains a type 1 thrombospondin domain, which is found in a number of proteins involved in the complement pathway, as well as in extracellular matrix proteins. Alternatively spliced transcript variants encoding different isoforms have been observed for this gene. [provided by RefSeq, Jan 2009]