

Product datasheet for RC212825L2

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc.

LGR5 (NM_003667) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: LGR5 (NM_003667) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: LGR5

Synonyms: FEX; GPR49; GPR67; GRP49; HG38

Mammalian Cell None

Selection:

Vector:pLenti-C-mGFP (PS100071)E. coli Selection:Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC212825).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_003667

ORF Size: 2721 bp





LGR5 (NM_003667) Human Tagged Lenti ORF Clone - RC212825L2

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 003667.2</u>

 RefSeq Size:
 2880 bp

 RefSeq ORF:
 2724 bp

 Locus ID:
 8549

 UniProt ID:
 075473

 UniProt ID:
 O75473

 Cytogenetics:
 12q21.1

Domains: 7tm 1, LRRNT, LRR, LRR TYP, LRR PS

Protein Families: Druggable Genome, GPCR, Transmembrane

MW: 99.8 kDa


Gene Summary: The protein encoded by this gene is a leucine-rich repeat-containing receptor (LGR) and

member of the G protein-coupled, 7-transmembrane receptor (GPCR) superfamily. The encoded protein is a receptor for R-spondins and is involved in the canonical Wnt signaling pathway. This protein plays a role in the formation and maintenance of adult intestinal stem cells during postembryonic development. Several transcript variants encoding different

isoforms have been found for this gene. [provided by RefSeq, Sep 2015]



Product images:



Circular map for RC212825L2