

Product datasheet for RC212874L4V

MAGEB10 (NM_182506) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MAGEB10 (NM_182506) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MAGEB10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_182506
ORF Size:	1041 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212874).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_182506.2 , NP_872312.2
RefSeq Size:	1953 bp
RefSeq ORF:	1044 bp
Locus ID:	139422
UniProt ID:	Q96LZ2
Cytogenetics:	Xp21.3
MW:	38.9 kDa
Gene Summary:	This gene encodes a member of the B subfamily of the melanoma associated antigen protein family. The encoded protein is specifically expressed in testis and tumor cells. [provided by RefSeq, Apr 2010]


[View online »](#)