

Product datasheet for **RC212876L4V**

FMNL3 (NM_198900) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FMNL3 (NM_198900) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FMNL3
Synonyms:	FHOD3; FRL2; WBP-3; WBP3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_198900
ORF Size:	2928 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212876).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_198900.2
RefSeq Size:	11039 bp
RefSeq ORF:	2931 bp
Locus ID:	91010
UniProt ID:	Q8IVF7
Cytogenetics:	12q13.12
MW:	111.1 kDa


[View online »](#)

Gene Summary:

The protein encoded by this gene contains a formin homology 2 domain and has high sequence identity to the mouse Wbp3 protein. Two alternative transcripts encoding different isoforms have been described. [provided by RefSeq, Jul 2008]