

## Product datasheet for **RC212885L3V**

### **KIR5.1 (KCNJ16) (NM\_170741) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	KIR5.1 (KCNJ16) (NM_170741) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KIR5.1
Synonyms:	BIR9; KIR5.1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_170741
ORF Size:	1254 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212885).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_170741.1</a> , <a href="#">NP_733937.1</a>
RefSeq Size:	4000 bp
RefSeq ORF:	1257 bp
Locus ID:	3773
UniProt ID:	<a href="#">Q9NPI9</a>
Cytogenetics:	17q24.3
Protein Families:	Druggable Genome, Ion Channels: Potassium, Transmembrane
MW:	48 kDa



[View online »](#)

**Gene Summary:**

Potassium channels are present in most mammalian cells, where they participate in a wide range of physiologic responses. The protein encoded by this gene is an integral membrane protein and inward-rectifier type potassium channel. The encoded protein, which tends to allow potassium to flow into rather than out of a cell, can form heterodimers with two other inward-rectifier type potassium channels. It may function in fluid and pH balance regulation. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Apr 2014]