

Product datasheet for RC212903L2

C14orf43 (ELMSAN1) (NM_194278) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	C14orf43 (ELMSAN1) (NM_194278) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	C14orf43
Synonyms:	C14orf43; C14orf117; c14_5541; ELMSAN1; LSR68
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212903).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.

ACCN:	NM_194278
ORF Size:	3135 bp



[View online »](#)

