

Product datasheet for **RC212939A1V**

Human LRRN4CL (NM_203422) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human LRRN4CL (NM_203422) AAV Particle
Tag: Myc-DDK
Symbol: LRRN4CL
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC212939 representing NM_203422
Red=Cloning site Blue=ORF Green=Tags(s)

CTATAGGGCGGCCGGAATTTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCCGCGCGC
GCC

ATGCTGGGCTCTCCCTGCCTTCTGTGGCTCCTGGCCGTGACCTTCTTGGTTCCAGAGCTCAGCCCTTGG
CCCCTCAAGACTTTGAAGAAGAGGAGGCAGATGAGACTGAGACGGCGTGGCCGCCTTTGCCGGCTGTCCC
CTGCGACTACGACCACTGCCGACACCTGCAGGTGCCCTGCAAGGAGCTACAGAGGGTCGGCCCGCGGCC
TGCCTGTGCCAGGACTCTCAGCCCCGCCAGCCGCCGACCCGCCGCGCATGGGAGAAGTGCGCATTG
CGGCCGAAGAGGGCCGCGCAGTGGTCCACTGGTGTGCCCTTCTCCCCGGTCTCCACTACTGGCTGCT
GCTTTGGGACGCGCAGGAGCTGCGCAGAAGGGGCCCGCTGAACGCTACGGTCCGAGAGCCGAAGT
AAGGGGCTGAAGCCAGGGGGCATTATGTCTGTTGCGTAGTGGCCGCTAACGAGGCCGGGCAAGCCGCG
TGCCCCAGGCTGGAGGAGAGGGCCTCGAGGGGGCCGACATCCCTGCCTTCGGGCCTTGACGCCGCTTGC
GGTGCCGCCCAACCCCGCACTCTGGTCCACGCGGCCGTGGGGTGGGCACGGCCCTGGCCCTGCTAAGC
TGTGCCGCCCTGGTGTGGCACTTCTGCCTGCGCATCGCTGGGGTGGCCGCGCCGAGCCGCCGCCGAG
CCGAGGGGCGCTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC212939 representing NM_203422
Red=Cloning site Green=Tags(s)

MLGSPCLLWLLAVTFLVPRAQPLAPQDFEEEEADETETAWPPLPAVPCDYDHCRLQVPCKELQRVGPAA
CLCPGLSSPAQPPDPPRMGEVRIAAEEGRAVVHWCAPFSPVLHYWLLLDWGSEAAQKGPPLNATVRRREL
KGLKPGGIYVVCVVAANEAGASRPQAGGEGLEGADIPAFGPCSRLAVPPNPRTLVAHAVGVGTALALLS
CAALVWHFCLDRWGCPRRAAARAAGAL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_203422
ORF Size:	714 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_203422.1</u>
RefSeq Size:	2429 bp
RefSeq ORF:	717 bp
Locus ID:	221091
UniProt ID:	<u>Q8ND94</u>
Cytogenetics:	11q12.3
MW:	25.1 kDa