

Product datasheet for **RC212943**

C10orf113 (NM_001010896) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C10orf113 (NM_001010896) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: C10orf113
Synonyms: bA165O3.1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC212943 representing NM_001010896
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTAAAAGTGAAAGAAGAATCTACTCTATGGAAAGCATGGCACCTGAAATTTCTGAAGATATTGGCT
GTCCTCTAGCTTATATGAGAGAGAGTGTGGTCTTACTAATCAACCAGTCATTTTTTTCTGTGTGGC
TGAAATGTACATTCCAGACATGAACAGGCCTGGTAAGGCAGAACGGAGAGCCGAGAAGGGCTGAGGCT
CCCAGGAAGCTGGTTGGGAGACCAGAAGGCTGCAAAGCTGGCCAACGCTTGCCAGGTGGCTGCGGAGAG
TCAGCTTCAGTCTCTACAAGCCGCCATCCAGGCAGCTCCTGAACCAGATCCAGGCAACTGGCCAAGTC
CCCAAGGCTGTCTCTGGGACCAGAAGGTATAACCAAAGGGAAGCATGGATTTAAATTCAGGCATAAAA
GAAAAGTTCACGTATCTAAGAAGGTGCTGAAAATGACTTTTTTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC212943 representing NM_001010896
Red=Cloning site Green=Tags(s)

MAKSERRIYSMESMAPEI SEDIGCLAYMRESV CALLINQSFFSCVAEMYIPDMNRPKAERRAEKGLRP
PRKLVGETRRLQSWPTLARWLRVRSFSLYKPIQAAPEPDGNWAKSPRLSLGPEGITKGKHKFKFQGIK
EKFNVS KVKLKM TFL

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV


Chromatograms: https://cdn.origene.com/chromatograms/mk8023_b02.zip



[View online >](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001010896

ORF Size: 465 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001010896.2](#), [NP_001010896.2](#)

RefSeq Size: 518 bp


RefSeq ORF: 468 bp

Locus ID: 387638

Cytogenetics: 10p12.31

MW: 17.5 kDa

Product images:



Circular map for RC212943